



## **Advantages**

Weighs Infants and Toddlers

Dual-purpose scale can be used with the infant weighing tray or toddler step-on platform.

Weight Hold

Weight value is held on the display so caregiver can attend to patient before recording the weight.

Imperial and Metric Measurement
Offers users the opportunity to weigh
patients in whichever unit is required.

**▼** Built-In Connection Port for eBABYHR

Digital baby height rod simply snaps onto the scale for a comprehensive weighing and measuring device.

Removable Tray

The removable tray makes it easy to clean the device.

**▼** Tare Function

Conveniently subtracts the weight of a blanket or diaper from the patient's weight.

**Easy to Use** 

Lightweight and portable with simple functionality and easy to read digital display.



## $\boldsymbol{\omega}$ S S





## **Specifications**

- Capacity: 44 lb / 20 kg
- Resolution: 0.5 oz / 10 g
- Tray Size (w x d x h):
   23 ½" x 8" x 2"
   597 mm x 203 mm x 51 mm
- Platform Size: (w x d x h):
   12 ½" x 8 ¾" x 1 ¼"
   318 mm x 219 mm x 32 mm
- Product Footprint Tray (w x d x h):
   23 ½" x 13 ¼" x 3 5%"
   597 mm x 337 mm x 92 mm
- Product Footprint Platform: (w x d x h):
   13 ½" x 11 ½" x 1 ¾"
   343 mm x 292 mm x 44 mm
- Display: ½" / 13 mm digit height LCD
- Power Source: 9V battery included
- Functions: Auto Off, Hold, Tare, Unit
- 2 Year Limited Warranty

Model	Description
SS-550KL	ScaleSurance 2 year extended warranty*
eBABYHR	Digital baby height rod for 550KL
	(must purchase with 550EHRBRACKET to mount on scale)
550EHRBRACKET	550KL mounting bracket for eBABYHR
	(must purchase with eBABYHR to mount on scale)
553CASE	Optional carrying case for 550KL
2210CART	Optional rolling cart for 550KL



Health o meter® Professional Scales 9500 West 55th Street USA McCook, IL 60525-7110 Phone: 1-800-815-6615

Fax: 1-800-654-7330

\* Not available in all countries